



Phonak Virto™ 1.

Product information

Phonak Virto Infinio is the comprehensive ITE portfolio for mild to profound hearing losses.

Powered by the innovative technology of the Phonak Infinio platform, Virto I hearing aids are the solution for those clients that prefer a discreet, custom-made solution.





Product description



Virto I-R
190/170/150/130



Virto I-10 NW O
190/170/150/130



Virto I-Titanium
190/170

Power level

	Virto I-R	Virto I-10 NW O	Virto I-Titanium
M	•	•	•
P	•	•	•
SP	•	•	•
UP ¹	•		

Hardware configuration / options / technology

	Virto I-R	Virto I-10 NW O	Virto I-Titanium
Battery size	built-in rechargeable	10	10
Charging times ³	0% to 100%	3h	
	0% to 80%	1.5h	
	30 min. fitting	15 min.	
Bluetooth®	•		
Shell material	Acrylic	Acrylic	Titanium*
IP rating	IP68 ²	IP68 ²	IP68 ²
Microphone technology	Directional	Omni	Omni
Telecoil		Optional	Optional
Push button (Program / Volume)	•	Optional	Optional
Volume control (wheel)		Optional	
EasyPhone (for acoustic phone)		•	Optional
MiniControl (for program or volume change)		•	•

Colors

Faceplate colors



26
Pink



14⁵
Tan



22
Cocoa



28⁵
Brown



06
Black

Titanium shell color



Titanium
gray

Acrylic shell colors



26
Pink



14
Tan



22
Cocoa



28
Brown



06
Black



R
Red Transparent



B
Blue Transparent



21
Transparent

¹ Available in specific countries only

² Actual charging time depends on the remaining battery life, but it will not exceed 3 hours.

³ Available in two different surface roughnesses: standard or extra retention (ERS)

⁴ IP68 indicates that the hearing aid is water and dust resistant. It survived continuous immersion in 1 meter of fresh water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard.

⁵ Virto I-R and Virto I-10 NW O only

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Sonova AG is under license.



Performance levels

AutoSense OS™ 6.0	Premium	Advanced	Standard	Essential
	I90	I70	I50	I30
Speech in Loud Noise ¹	•			
Music	•	•		
Speech in Car	•	•	•	
Comfort in Noise	•	•	•	
Speech in Noise	•	•	•	•
Calm Situation	•	•	•	•
Media speech + mic ¹	•	•	•	•
Media music + mic ¹	•	•	•	•
Additional programs				
Max. additional programs	4	4	4	4
Speech in Loud Noise ¹	•	•		
Speech in 360°	•	•		
Speech in car	•	•	•	
Speech in Noise	•	•	•	•
Calm Situation	•	•	•	•
Comfort in Noise	•	•	•	•
Music	•	•	•	•
Acoustic phone	•	•	•	•
Custom program	•	•	•	•
Phone via T-coil ± mic ²	•	•	•	•
Public T-coil ± mic ²	•	•	•	•
Streaming programs				
RogerDirect™ + mic ¹	•	•	•	•
PartnerMic + mic ¹	•	•	•	•
Phone call + mic ¹	•	•	•	•

Features	Premium	Advanced	Standard	Essential
	I90	I70	I50	I30
SpeechSensor ¹	•			
Speech Enhancer	•			
StereoZoom 2.0 ¹	•	•		
Dynamic Noise Cancellation ¹	•	•		
Motion Sensor Hearing ¹	•	•	•	
Roger compatibility	•	•	•	•
Tap control ¹				
- Siri/Google Assistant	•	•		
- Pause/resume streaming	•	•		
- Accept/end phone call	•	•		
SoundRelax	•	•	•	
DuoPhone ¹	•	•	•	
EasyPhone ²	•	•	•	•
UltraZoom ¹	•	•	•	•
Roger directional setting ¹	•	•	•	•
WindBlock	•	•	•	•
User Preference Tuning	•	•	•	•
Tinnitus Balance	•	•	•	•
Auto Acclimatization	•	•	•	•
Environmental balance	•	•	•	•
SoundRecover2	•	•	•	•
WhistleBlock	•	•	•	•
NoiseBlock	•	•	•	•
QuickSync ¹	•	•	•	•
AOV	•	•	•	•
Fine tuning channels	20	20	16	12





¹ Available for Phonak Virto I-R only
² Available for Phonak Virto I-10 NW 0 and Virto I-Titanium only

Fitting

	Virto I-R	Virto I-10 NW 0 / Virto I-Titanium
Phonak Target	10.0 or higher	10.0 or higher
Noahlink Wireless, Noahlink Wireless 2	•	
NOAHlink		•
HI-PRO, HI-PRO2		•



Acoustic coupling

		M	P	SP	UP ¹
Maximum Power Output (dB SPL)	2 cc coupler	109	115	119	127
Maximum gain (dB)	2 cc coupler	40	50	60	70
Frequency range (Hz)	2 cc coupler	<100 - 7000	<100 - 6700	<100 - 7000	<100 - 5200
Working current (mA)	Virto I-10 NW O, Virto I-Titanium	1.1	1.2	1.1	
Fitting range					
		For mild to moderate hearing loss, all audiometric configurations	For mild to moderately-severe hearing loss, all audiometric configurations	For moderate to severe hearing loss, all audiometric configurations	For moderate to profound hearing loss, all audiometric configurations

¹ Available for Virto I-R and in specific countries only

Shell styles



Virto I-R



Virto I-10 NW O



Virto I-Titanium

	Virto I-R	Virto I-10 NW O	Virto I-Titanium
IIC		•	•
CIC		•	•
MC		•	
ITC	•	•	
HS	•	•	
FS	•	•	

For more information please visit www.phonak.com/virto-infinio
 For information on Roger microphones please visit www.phonak.com/roger
 For more information on myPhonak please visit www.phonak.com/apps





Phonak Virto I-R charging options

Phonak Charger Virto I

Phonak Virto I-R has all the Infinio features you love as well as the benefits of rechargeability.

Clients can enjoy a full day of hearing including streaming*.



* Expected battery life of 21 hours, including 4 hours of BT classic streaming.

Phonak Virto I-R Bluetooth® version

Direct calls/audio streaming

Bluetooth 5.2 wireless technology and most older Bluetooth phones

Wireless communication portfolio

	Virto I-R 190/170/150/130	Virto I-10 NW 0 190/170/150/130	Virto I-Titanium 190/170
Wireless accessory			
TV Connector	•		
Phonak PartnerMic	•		
Phonak RemoteControl	•		
Roger			
RogerDirect	•		
Roger Neckloop		•*	•*
Apps			
myPhonak	•		

* Available if telecoil is installed





life is on


Since 1947, Phonak is dedicated to preserving life quality by opening new acoustic worlds. Already back then, in the days of our foundation, our company was driven by a formative conviction: We believe that well-hearing equates to well-being and thus is essential for living life to the fullest. In fact, the sense of hearing is directly linked to social, emotional, cognitive and physical well-being. Today as in future, we thrive to offer the broadest portfolio of innovative hearing solutions. And, together with our hearing care professionals, we keep on focusing on what matters most: improving speech understanding, changing people's lives and having a positive effect on society as a whole.

phonak.com



027-0710-02/M1.00/2024-08/na © 2024 Sonova AG All rights reserved



 Sonova AG · Laubisrütistrasse 28
CH-8712 Stäfa · Switzerland

